

CR39 FLAT TOP



VIDCON

CR39 FLAT TOP

VIDEON

Material refractive index=1.499

Nominal Curve	Dia	Radius of Curvature	True Curve	Sag	Back Curve	Seg inset	Seg Drop	Add Range	C.T.	E.T.	Rx Range Recommended
50	70/28mm	981,48	0,54	0,32	5,70	5	5	0.50-4.00	10,00	16,00	-5.25D/-10.00D
200	70/28mm	262,38	2,02	1,19	5,70	5	5	0.50-6.00	11,60	16,00	-2.00D/-5.00D
400	70/28mm	135,20	3,92	2,33	5,70	5	5	0.50-6.00	6,80	9,00	-1.00D/-3.00D
500	70/28mm	104,74	5,06	3,03	5,70	5	5	0.75-4.00	8,40	9,00	+1.00D/-2.00D
600	70/28mm	87,46	6,06	3,65	5,70	5	5	0.50-6.00	9,50	9,00	+2.00D/-2.00D
700	70/28mm	74,96	7,07	4,29	5,70	5	5	0.50-4.00	10,50	9,00	+0.00D/+4.00D
800	70/28mm	65,84	8,05	4,93	5,70	5	5	0.50-6.00	10,00	7,00	+2.00D/+5.00D
1000	70/28mm	52,53	10,09	6,33	5,70	5	5	0.50-4.00	13,00	7,00	+5.25D/+7.00D

1.499 S/F FLAT TOP BIFOCAL LENS (70/28mm)